



***Filling in the 'Green Gap' LED challenge for  
microscopy and  
medical illumination. Providing mercury-  
free sustainable illumination***

Dr. Kavita Aswani, PhD

Senior Applications Scientist

IS2L Denver October 2019


# LED vs. Lamp



LED	Lamp
No flicker	Flicker
Instant ON/OFF	Warm-up, cool-down time
Slower intensity loss over time	30% intensity loss in first 400 hours
More stable	Less so
Repeatable output power	
Long lifetime	Lower lifetime, need bulb stock
Less phototoxic to cells	
69% less power consumption	
Zero Mercury	Larger green footprint
No consumables	Bulbs and Light guides
Selectable wavelengths	Reliant on a filter wheel



# LED vs. Lamp intensity over time



## Cost of Ownership (per 20,000 hours of "ON time")

	HBO	X-Cite 120Q	X-Cite 120LED <i>mini</i>	X-Cite XYLIS
Replacement Lamps	100	10	-	-
Mercury Content	1 100 mg	200 mg	-	-
Lamp Costs	\$15,000	\$6,250	-	-
Replacement Light Guides	-	5	-	2
Light Guide Costs	-	\$1,975	-	\$790
Bulb Disposal (\$5/bulb) <sup>1</sup>	\$500	\$50	-	-
Maintenance Costs (bulb, \$20/hr) <sup>2</sup>	\$1,000	\$17	-	-
<b>TOTAL</b>	<b>\$16,500</b>	<b>\$8,292</b>	<b>\$0</b>	<b>\$790</b>
<b>Hourly Cost<sup>3</sup></b>	<b>\$0.82</b>	<b>\$0.41</b>	<b>\$0</b>	<b>\$0.04</b>
<b>Annual Cost<sup>4</sup></b>	<b>\$1,650</b>	<b>\$829</b>	<b>\$0</b>	<b>\$10</b>

### Notes:


1. Mercury-Free Microscopy white paper [www.mygreenlab.org](http://www.mygreenlab.org).
2. Assumes 30 min to change/align HBO lamp, 5 min for X-Cite 120Q.
3. Assumes 8 hour day, 4x15 min. imaging sessions. Arc lamps left on for the day and LEDs on continuously during each session.
4. Assumes a 5 day week x 50 weeks.




# What's the issue?



- Fluors developed using lamp spectrum
  - No green gap in lamps
- LEDs sufferED from the Green Gap
- LEDs WERE not bright enough for FISH




# What is the Green Gap?




- Lack of Semiconductor material that emits light between 540-590nm
- Issue for fluors – green/yellow excitation, red emission
- Lamps and Lasers used traditionally

By User:S-kei - File:PnJunction-LED-E.PNG, CC BY-SA 2.5, <https://commons.wikimedia.org/w/index.php?curid=14985902>


# Lamp Output in the Green Gap




# Then came LED











# LaserLED Hybrid Drive™




## Laser phosphor conversion

Image courtesy of [www.lightemittingdiodes.com](http://www.lightemittingdiodes.com)

# Solving the Green Gap challenge



# Spectral Comparison – LED vs. Lamp




## Notes:

- X-Cite XYLIS Prototype Data - Beta version, August 2017
- For LLG coupled systems, measurements were obtained with 3mm x 1.5m Liquid Light Guide
- For direct coupled systems, measurements were scaled down for realistic comparison to LLG coupled systems
- In all cases, light coupling efficiency depends on the microscope optics
- Data represents typical output levels. Output will vary between individual units

# Power Comparisons



Power at Specimen Plane  
Olympus BX50, 10X objective




## Notes:

1. X-Cite XYLIS Prototype – Beta version, August 2017
2. Measurements taken with XR2100 & XP750
3. Light sources coupled to microscope with 810-00032X adaptor
4. Excitation Filters: 560/25, 650/13, 710/75

# Power Comparisons



Power at Specimen Plane  
Zeiss AxioObserver, 20X objective



## Notes:

1. X-Cite XYLIS Prototype – Beta version, July 2017
2. Measurements taken with XR2100 & XP750
3. Light sources coupled to microscope with 810-00022X adaptor
4. Excitation Filters: 375/25, 470/40, 545/25, 560/40, 620/60, 710/75

# Why LED?



- **Less down time – no bulbs to replace**
- **NOW enough power for FISH**
- **Uniform and stable**
- **Repeatable exposure**
- **Faster throughput with switchable LEDs**
- **Less waste/toxic disposal**
- **Less energy consumption**
- **Less running costs**



# Rejuvenated LED applications



## Analytical Solutions



### APPLICATIONS

- Fluorescence microscopy/microfluidics
- Fluorescence slide/plate scanners
- Microscopy, DMD

## Medical Solutions




### APPLICATIONS

- Endoscopy
- Surgical Visualisation

# Vitae vIR spectra



	Radiant Flux (W)	CCT
300W Xenon	4.8	5110K
vIR Boost	4.25	5200K
vIR	2.42	5100K
White LED	1.34	5640K



# Conclusion:



- LEDs no longer have to suffer from Green Gap
- LaserLED Hybrid drive provide the closest spectrum to lamp output
- Rest easy with no lamp replacements, cost and disposal
- Keep your fluorophores and cells happy 😊



**Thank you for your attention.**